

Product Datasheet

ALB Antibody (Biotin) (orb243722)

Catalog Number	orb243722
Category	Antibodies
Description	This ALB Antibody, Biotin conjugated is an Biotin polyclonal product. It targets IgG using a Serum Albumin purified from Rabbit serum (Native Protein) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified.
Target	ALB
Clonality	Polyclonal
Species/Host	Guinea pig
Isotype	IgG
Conjugation	Biotin
Reactivity	Rabbit
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Serum Albumin purified from Rabbit serum (Native Protein)
UniProt ID	P49065
Tested applications	ELISA
Dilution range	IHC-P:1:20-200, IF:1:50-200

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)