

Product Datasheet

AMPK beta 1/PRKAB1 Rabbit Polyclonal Antibody (PerCP) (orb2461695)

Catalog Number	orb2461695
Description	AMPK beta 1/PRKAB1 Rabbit Polyclonal Antibody (PerCP)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PerCP
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human AMPK beta 1 (2-100/270aa)
Target	PRKAB1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	ICC/IF=1:50, IF=1:100-500, Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)