

Product Datasheet

Yeast PTPRZ1 protein (orb246780)

Catalog Number	orb246780
Category	Proteins
Description	This Yeast PTPRZ1 protein spans the amino acid sequence from region 36-300aa. Purity: Greater than 90% as determined by SDS-PAGE.
Tag	N-terminal 6xHis-tagged
Form/Appearance	Liquid or Lyophilized powder
Buffer/Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Purity	Greater than 90% as determined by SDS-PAGE.
Protein Sequence	IGWSYTGALNQKNWGKKYPTCNSPKQSPINIDEDLTQVNVNLKKLKFQGWDKTSLENTFIH NTGKTVEINLTNDYRVSGGVSEMVFKASKITFWHGKCNMSSDGSEHSLEGQKFPLEMQIYC FDADRFSSFEEAVKGKGKLRALSILFEVGTEENLDFKAIIDGVESVSRFGKQAALDPFILLNLL PNSTDKYYIYNGSLTSPPCDTVDWIVFKDTSISESQLAVFCEVLTMQQSGYVMLMDYLQ NNFREQQYKFSRQVFSSY
Protein Length	Partial
UniProt ID	P23471
MW	32.1 kDa
Application notes	This is His-tag protein
Source	Yeast
Biological Origin	Homo sapiens (Human)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expression Region


36-300aa

Storage


The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note

For research use only



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of Yeast host and target protein, the LC-MS/MS Analysis result could indicate that this peptide derived from Yeast-expressed Homo sapiens (Human) PTPRZ1.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Based on the SEQUEST from database of Yeast host and target protein, the LC-MS/MS Analysis result could indicate that this peptide derived from Yeast-expressed Homo sapiens (Human) PTPRZ1.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558