

Product Datasheet

NKp46 Rabbit Polyclonal Antibody (FITC) (orb247242)

Catalog Number	orb247242
Category	Antibodies
Description	NKp46 Rabbit Polyclonal Antibody (FITC) is an FITC conjugated polyclonal antibody that specifically recognizes NKp46. It is supplied in liquid form and exhibits reactivity with Human, Mouse and Rat samples. The suitable applications for this antibody are ICC and IF.
Target	NKp46 (FITC)
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Human, Mouse, Rat
Form/Appearance	0.5 mg/ml in buffer 1*PBS; or 0.2 mg/ml in PBS based buffer containing 30% glycerol and 0.1% sodium
Concentration	batch dependent
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
UniProt ID	Q8C567, Q9Z0H5, O76036
Tested applications	ICC, IF
Dilution range	IF/ICC: 1:100

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Entrez	9437
NCBI	NP_001138929.2, NP_001138929.21, NM_001145457.2, NM_001145457.21
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of human hela cells tissue using NKp46 antibody (dilution of primary antibody - 1:100)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)