

Product Datasheet

Phospho-MST4 + MST3 + STK25 (Thr178 + Thr190 + Thr174) Rabbit Polyclonal Antibody (PerCP-Cy5.5) (orb2552202)

Catalog Number	orb2552202
Category	Antibodies
Description	Phospho-MST4 + MST3 + STK25 (Thr178 + Thr190 + Thr174) Rabbit Polyclonal Antibody (PerCP-Cy5.5) is a PerCP/Cy5.5 conjugated antibody targeting STK25. This antibody is suitable for IF. It exhibits reactivity with Human, Mouse, Rat samples.
Target	STK25
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PerCP/Cy5.5
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Equine, Sheep
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Immunogen	KLH conjugated synthesised phosphopeptide derived from human MST4/MST3/STK25 around the phosphorylation site of Thr178/Thr190/Thr174 RN(p-T)FV
Tested applications	IF
Dilution range	IF=1:100-500
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)