

## Product Datasheet

### Rabbit IgG, BF555 conjugated (orb2560635)

|                             |  |
|-----------------------------|--|
| <b>Catalog Number</b>       | orb2560635   |
| <b>Category</b>             | Proteins   |
| <b>Description</b>          | This product is Rabbit IgG, BF555 conjugated, also known as Rabbit Immunoglobulin G/BF555 Isotype Control. For optimal stability, store at -20°C protected from light and avoid repeated freeze-thaw cycles. |
| <b>Form/Appearance</b>      | Liquid   |
| <b>Buffer/Preservatives</b> | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol.   |
| <b>Storage</b>              | Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  |
| <b>Note</b>                 | For research use only  |
| <b>Expiration Date</b>      | 6 months from date of receipt.   |

---

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)