

Product Datasheet

Bovine Pepsinogen A (PGA) ELISA Kit (orb2561063)

Catalog Number	orb2561063
Category	Assays and Kits
Description	The Bovine Pepsinogen A (PGA) ELISA Kit features a dynamic range of 1.56-100ng/mL and is suitable for Serum, plasma or other biological fluids. This product offers a sensitivity of 0.58ng/mL.
Reactivity	Bovine
Range	1.56-100ng/mL
Sensitivity	0.58ng/mL
Sample Types	Serum, plasma or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)