

## Product Datasheet

### RANTES Recombinant Rabbit Monoclonal Antibody pair (capture) (orb2563366)

|                     |                                                                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number      | orb2563366                                                                                                                              |
| Description         | RANTES Recombinant Rabbit Monoclonal Antibody pair (capture)                                                                            |
| Species/Host        | Rabbit                                                                                                                                  |
| Reactivity          | Human                                                                                                                                   |
| Conjugation         | Unconjugated                                                                                                                            |
| Tested Applications | ELISA                                                                                                                                   |
| Immunogen           | Recombinant protein                                                                                                                     |
| Target              | CCL5                                                                                                                                    |
| Preservatives       | 0.01M PBS, pH 7.2.                                                                                                                      |
| Form/Appearance     | Liquid                                                                                                                                  |
| Storage             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note                | For research use only                                                                                                                   |
| Isotype             | IgG                                                                                                                                     |
| Clonality           | Recombinant                                                                                                                             |
| Antibody Type       | Recombinant Antibody                                                                                                                    |
| Expiration Date     | 12 months from date of receipt.                                                                                                         |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)