

Product Datasheet

Rabbit Anti-Goat IgG H&L, BF555 conjugated (orb2563570)

Catalog Number	orb2563570
Category	Antibodies
Description	Rabbit Anti-Goat IgG H&L, BF555 conjugated
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	BF555
Reactivity	Goat
Concentration	2.0 mg/ml
Buffer/Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Purification	Affinity purified by Protein A
Immunogen	Native Goat IgG
Tested applications	FC, ICC, IF
Dilution range	IF=1:100-1000, ICC/IF=1:100-1000, Flow-Cyt=1:100-1000
Antibody Type	Secondary Antibody
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558