

Product Datasheet

Anti-V5 epitope tag [SV5-P-K] (orb256435)

Catalog Number	orb256435
Category	Antibodies
Description	Human monoclonal antibody to V5 tag
Target	V5 epitope tag
Clonality	Monoclonal
Species/Host	Mouse
Isotype	Mouse Fab fragment
Conjugation	Unconjugated
Reactivity	Virus
Concentration	1 mg/ml
Buffer/Preservatives	PBS with 0.02% Proclin 300.
Purity	Purified
Immunogen	Pk (V5) Epitope Tag (GKPIPPLLGLDST).
Tested applications	ELISA, IF, IHC-Fr, IP, RIA, WB
Specificity	Recognises a 14 amino acid sequence (GKPIPPLLGLDST) derived from the P and V proteins of the paramyxovirus simian virus 5.
Clone Number	SV5-P-K
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.



ELISA of anti-epitope tag antibodies on multi-tag protein Plate coated with 5 ug/ml of Multitag protein. Primary antibody orb256432 loaded at 1000 ng/ml in first well and an 11 point serial dilution performed in duplicate for each antibody with last well left blank. Signal detection was performed using HRP-conjugated goat anti-mouse.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558