


Product Datasheet

Human MIP-3 α /CCL20 (Macrophage Inflammatory Protein 3 Alpha) Quick ELISA Kit (orb2568539)

Description	Human MIP-3 α /CCL20 (Macrophage Inflammatory Protein 3 Alpha) Quick ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	MIP-3 α
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Assay Time	120 minutes
Uniprot ID	P78556
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)