


Product Datasheet

Human CTSD (Cathepsin D) ELISA Kit (orb2568689)

Catalog Number	orb2568689
Category	Assays and Kits
Description	Human CTSD (Cathepsin D) ELISA Kit
Target	CTSD
Reactivity	Human
UniProt ID	P07339
Assay Type	Sandwich
Assay Time	4 hours
Range	0.313-20ng/ml
Sensitivity	0.188ng/ml
Storage	2 - 8 °C for 12 months
Note	For research use only



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558