

## Product Datasheet

### JUNB (phospho-S259) peptide (orb256931)

<b>Catalog Number</b>	orb256931
<b>Category</b>	Proteins
<b>Description</b>	The peptide is used to block Anti-JUNB (pS259) Antibody (orb256641) reactivity.
<b>Target</b>	JUNB
<b>Reactivity</b>	Bovine, Human, Mouse, Rat, Sheep
<b>Form/Appearance</b>	Lyophilized powder
<b>Purity</b>	>85%
<b>UniProt ID</b>	<b>P17275, P24898, P09450</b>
<b>Tested applications</b>	Blocking
<b>Application notes</b>	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Source</b>	Synthetic
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Entrez</b>	<b>3726, 16477, 24517</b>
<b>Expiration Date</b>	12 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)