

Product Datasheet

LOX-1 Rabbit Polyclonal Antibody (FITC) (orb257059)

Catalog Number	orb257059
Category	Antibodies
Description	LOX-1 Rabbit Polyclonal Antibody (FITC)
Target	LOX-1 (FITC)
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Form/Appearance	0.5 mg/ml in buffer 1*PBS; or 0.2 mg/ml in PBS based buffer containing 30% glycerol and 0.1% sodium
Concentration	batch dependent
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Immunogen	Synthetic peptide (please contact us for details)
UniProt ID	Q9EQ09, P78380, O70156
MW	31 kDa
Tested applications	ICC, IF

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Entrez	4973
NCBI	NP_001166103.1, NM_001172632.11, NM_001172632.1, NP_001166103.11
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)