

Product Datasheet

Anti-Cuscuta campestris CCAM_LOCUS22811 Antibody (orb2572088)

Description	Anti-Cuscuta campestris CCAM_LOCUS22811 Antibody. This antibody reacts with Cuscuta campestris.
Species/Host	Rabbit
Reactivity	Plant
Conjugation	APC
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Uniprot ID	A0A484LXA6
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)