

Product Datasheet

Anti-Cuscuta campestris CCAM_LOCUS22811 Antibody (orb2572092)

Catalog Number	orb2572092
Category	Antibodies
Description	Anti-Cuscuta campestris CCAM_LOCUS22811 Antibody. This antibody reacts with Cuscuta campestris.
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Rabbit IgG
Conjugation	Biotin
Reactivity	Plant
Form/Appearance	Liquid
Concentration	500 µg/ml
Purification	Immunogen affinity purified.
UniProt ID	A0A484LXA6
Antibody Type	Primary Antibody
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)