

Product Datasheet

Anti-Trypanosoma brucei brucei Tb927.3.4760 Antibody (orb2572742)

Description	Polyclonal Anti-Trypanosoma brucei brucei (strain 927/4 GUTat10.1) Tb927.3.4760 Antibody. This antibody reacts with Trypanosoma brucei brucei (strain 927/4 GUTat10.1).
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Cy3
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Uniprot ID	Q582Q9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)