

Product Datasheet

Human NKp46 Protein, Fc Tag (orb257704)

Catalog Number	orb257704
Category	Proteins
Description	Human NKp46, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Asn 254 (Accession # AAH64806).
Target	NKp46
Conjugation	Unconjugated
Tag	C-Fc
Reactivity	Human
Form/Appearance	Powder
Buffer/Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 10 mM NaCl, pH 7.5
Purity	95%
UniProt ID	AAH64806
MW	53.0 kDa
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 53.0 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Expression System	HEK293
Endotoxins	0.01 EU per µg
Expression Region	Gln 22 - Asn 254

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage	-20°C
Note	For research use only
NCBI	AAH64806
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human NKp46 protein



Immobilized Human NKp46, Fc Tag, premium grade (Cat. No. orb257704) at 1 µg/mL (100 µL/well) can bind Anti-NKP46 Antibody, Human IgG1 with a linear range of 0.2-4 ng/mL (QC tested).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558