

Product Datasheet

Human Semaphorin 4A Protein, His Tag (orb257835)

Catalog Number	orb257835
Category	Proteins
Description	Human Semaphorin 4A, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 33 - His 683 (Accession # AAH20974).
Target	Semaphorin 4A
Conjugation	Unconjugated
Tag	C-6×His
Reactivity	Human
Form/Appearance	Powder
Buffer/Preservatives	PBS, pH 7.4
Purity	95%
UniProt ID	AAH20974
MW	72.7 kDa
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 72.7 kDa. The protein migrates as 70-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Expression System	HEK293
Endotoxins	1.0 EU per µg

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expression Region


Gly 33 - His 683

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Note

For research use only



SDS-PAGE analysis of Human Semaphorin 4A protein

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.comPhone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)**Biorbyt LLC**

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.comPhone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)