

## Product Datasheet

### Anti-NOXA1 Antibody (orb2580880)

|                            |                                                                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-NOXA1 Antibody. Tested in WB, ELISA applications. This antibody reacts with Human.                                                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                                                                                                             |
| <b>Conjugation</b>         | APC                                                                                                                                                                                               |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                                                                         |
| <b>Immunogen</b>           | E.coli-derived human NOXA1 recombinant protein (Position: A38-Q434). Human NOXA1 shares 60.3% and 61% amino acid (aa) sequence identity with mouse and rat NOXA1, respectively.                   |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                       |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                         |
| <b>Storage</b>             | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing. |
| <b>Note</b>                | For research use only                                                                                                                                                                             |
| <b>Application notes</b>   | Western blot, 0.25-0.5 µg/ml, Human ELISA, 0.1-0.5 µg/ml, -. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                        |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                                                                                  |
| <b>MW</b>                  | 60 kDa                                                                                                                                                                                            |
| <b>Uniprot ID</b>          | <b>Q86UR1</b>                                                                                                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)