

## Product Datasheet

### MIP-3 Beta/CCL19/MIP Rabbit Polyclonal Antibody (FITC) (orb2585187)

|                      |                                                                                                                     |
|----------------------|---------------------------------------------------------------------------------------------------------------------|
| Catalog Number       | orb2585187                                                                                                          |
| Category             | Antibodies                                                                                                          |
| Description          | Anti-MIP-3 Beta/CCL19 Antibody. Tested in ELISA, IHC, Flow Cytometry applications. This antibody reacts with Human. |
| Target               | C-C motif chemokine 19                                                                                              |
| Clonality            | Polyclonal                                                                                                          |
| Species/Host         | Rabbit                                                                                                              |
| Isotype              | IgG                                                                                                                 |
| Conjugation          | FITC                                                                                                                |
| Reactivity           | Human                                                                                                               |
| Form/Appearance      | Liquid                                                                                                              |
| Buffer/Preservatives | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .        |
| Purification         | Immunogen affinity purified.                                                                                        |
| Immunogen            | E.coli-derived human MIP-3 Beta/CCL19 recombinant protein (Position: G22-S98).                                      |
| UniProt ID           | <b>Q99731</b>                                                                                                       |
| MW                   | 52588 Da                                                                                                            |
| Dilution range       | Optimal dilutions should be determined by end users.                                                                |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                          |                                                                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Application notes</b> | Recommended applications are based on the parent unconjugated antibody. Customers may select suitable applications according to their experimental needs.                                         |
| <b>Cross Reactivity</b>  | No cross reactivity with other proteins.                                                                                                                                                          |
| <b>Antibody Type</b>     | Primary Antibody                                                                                                                                                                                  |
| <b>Storage</b>           | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing. |
| <b>Note</b>              | For research use only                                                                                                                                                                             |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)