

Product Datasheet

Anti-F4/80/Adgre1 Antibody (orb2586184)

Description	Anti-F4/80/Adgre1 Antibody. Tested in IF, IHC applications. This antibody reacts with Mouse, Rat.
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	iFluor647
Tested Applications	IF, IHC
Immunogen	E.coli-derived mouse F4/80/Adgre1 recombinant protein (Position: G27-S230).
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Mouse, Rat Immunofluorescence, 5 µg/ml, Mouse, Rat
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
MW	66 kDa
Uniprot ID	Q61549
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)