

## Product Datasheet

# Ubc13/UBE2N Rabbit Polyclonal Antibody (HRP) (orb2592844)

|                             |                                                                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2592844                                                                                                                         |
| <b>Category</b>             | Antibodies                                                                                                                         |
| <b>Description</b>          | Anti-Ubc13/UBE2N Antibody. Tested in ELISA, Flow Cytometry, IF, ICC, WB applications. This antibody reacts with Human, Mouse, Rat. |
| <b>Target</b>               | Ubiquitin-conjugating enzyme E2 N                                                                                                  |
| <b>Clonality</b>            | Polyclonal                                                                                                                         |
| <b>Species/Host</b>         | Rabbit                                                                                                                             |
| <b>Isotype</b>              | Rabbit IgG                                                                                                                         |
| <b>Conjugation</b>          | HRP                                                                                                                                |
| <b>Reactivity</b>           | Human, Mouse, Rat                                                                                                                  |
| <b>Form/Appearance</b>      | Liquid                                                                                                                             |
| <b>Buffer/Preservatives</b> | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .                                                |
| <b>Purification</b>         | Immunogen affinity purified.                                                                                                       |
| <b>Immunogen</b>            | E.coli-derived human Ubc13/UBE2N recombinant protein (Position: Q13-I152).                                                         |
| <b>UniProt ID</b>           | <b>P61088</b>                                                                                                                      |
| <b>MW</b>                   | 38779 Da                                                                                                                           |
| <b>Tested applications</b>  | ELISA, IHC, WB                                                                                                                     |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                         |                                                                                                                                                                                                                                       |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Dilution range</b>   | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |
| <b>Cross Reactivity</b> | No cross-reactivity with other proteins.                                                                                                                                                                                              |
| <b>Antibody Type</b>    | Primary Antibody                                                                                                                                                                                                                      |
| <b>Storage</b>          | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.                                                                                                                                                      |
| <b>Note</b>             | For research use only                                                                                                                                                                                                                 |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                                                                                                                                                                                       |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)