

## Product Datasheet

### Anti-CSNK2A2 Antibody (orb2597468)

|                            |                                                                                                                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-CSNK2A2 Antibody. Tested in Flow Cytometry, IF, ICC, WB applications. This antibody reacts with Human, Mouse, Rat.                                                                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                                                                                 |
| <b>Conjugation</b>         | iFluor647                                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | FC, ICC, IF, WB                                                                                                                                                                                                                                   |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence at the C-terminus of human CSNK2A2, which shares 93.3% amino acid (aa) sequence identity with mouse CSNK2A2.                                                                                      |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                                       |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                                                                         |
| <b>Storage</b>             | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                                 |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                             |
| <b>Application notes</b>   | Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat<br>Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry (Fixed), 1-3 µg/1x10 <sup>6</sup> cells, Human. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                        |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                                                                                                                                  |
| <b>MW</b>                  | 39 kDa                                                                                                                                                                                                                                            |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

|                        |                                 |
|------------------------|---------------------------------|
| <b>Uniprot ID</b>      | <b>P19784</b>                   |
| <b>Expiration Date</b> | 12 months from date of receipt. |

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)