

## Product Datasheet

# Kinectin 1/KTN1 Rabbit Polyclonal Antibody (HRP) (orb2599844)

|                      |                                                                                                                                                     |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number       | orb2599844                                                                                                                                          |
| Category             | Antibodies                                                                                                                                          |
| Description          | Anti-Kinectin 1/KTN1 Antibody. Tested in ELISA, Flow Cytometry, IF, IHC, ICC, WB applications. This antibody reacts with Human, Mouse, Rat, Monkey. |
| Target               | Kinectin                                                                                                                                            |
| Clonality            | Polyclonal                                                                                                                                          |
| Species/Host         | Rabbit                                                                                                                                              |
| Isotype              | Rabbit IgG                                                                                                                                          |
| Conjugation          | HRP                                                                                                                                                 |
| Reactivity           | Human, Monkey, Mouse, Rat                                                                                                                           |
| Form/Appearance      | Liquid                                                                                                                                              |
| Buffer/Preservatives | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .                                                                 |
| Purification         | Immunogen affinity purified.                                                                                                                        |
| Immunogen            | E.coli-derived human Kinectin 1/KTN1 recombinant protein (Position: K283-K827).                                                                     |
| UniProt ID           | <b>Q86UP2</b>                                                                                                                                       |
| MW                   | 39125 Da                                                                                                                                            |

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                         |                                                                                                    |
|-------------------------|----------------------------------------------------------------------------------------------------|
| <b>Dilution range</b>   | The intended application should be selected according to the customer's experimental requirements. |
| <b>Cross Reactivity</b> | No cross-reactivity with other proteins.                                                           |
| <b>Antibody Type</b>    | Primary Antibody                                                                                   |
| <b>Storage</b>          | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.                   |
| <b>Note</b>             | For research use only                                                                              |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                                                    |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)