

## Product Datasheet

### Anti-Arp3/ACTR3 Antibody (orb2605931)

|                            |                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-Arp3/ACTR3 Antibody. Tested in ELISA, Flow Cytometry, WB applications. This antibody reacts with Human, Mouse, Rat.                                                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                                                   |
| <b>Conjugation</b>         | iFluor647                                                                                                                                                                                                           |
| <b>Tested Applications</b> | ELISA, FC, WB                                                                                                                                                                                                       |
| <b>Immunogen</b>           | E.coli-derived human Arp3/ACTR3 recombinant protein (Position: M1-S418).                                                                                                                                            |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                         |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                                           |
| <b>Storage</b>             | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.               |
| <b>Note</b>                | For research use only                                                                                                                                                                                               |
| <b>Application notes</b>   | Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat Flow Cytometry (Fixed), 1-3µg/1x10 <sup>6</sup> cells, Human, Mouse, Rat ELISA, 0.1-0.5µg/ml, -. Add 0.2ml of distilled water will yield a concentration of 500ug/ml |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                          |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                                                                                                    |
| <b>MW</b>                  | 47 kDa                                                                                                                                                                                                              |
| <b>Uniprot ID</b>          | <b>P61158</b>                                                                                                                                                                                                       |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Expiration Date**

12 months from date of receipt.

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)