

## Product Datasheet

### Prealbumin/TTR Rabbit Polyclonal Antibody (Biotin) (orb2619439)

|                      |                                                                                                                                                              |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number       | orb2619439                                                                                                                                                   |
| Category             | Antibodies                                                                                                                                                   |
| Description          | Prealbumin/TTR Rabbit Polyclonal Antibody (Biotin)                                                                                                           |
| Target               | Transthyretin                                                                                                                                                |
| Clonality            | Polyclonal                                                                                                                                                   |
| Species/Host         | Rabbit                                                                                                                                                       |
| Isotype              | Rabbit IgG                                                                                                                                                   |
| Conjugation          | Biotin                                                                                                                                                       |
| Reactivity           | Mouse, Rat                                                                                                                                                   |
| Form/Appearance      | Liquid                                                                                                                                                       |
| Buffer/Preservatives | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .                                                 |
| Purification         | Immunogen affinity purified.                                                                                                                                 |
| Immunogen            | E.coli-derived rat Prealbumin recombinant protein (Position: G21-N147). Rat Prealbumin shares 94.5% amino acid (aa) sequence identity with mouse Prealbumin. |
| UniProt ID           | <b>P02767</b>                                                                                                                                                |
| MW                   | 15720 Da                                                                                                                                                     |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                         |                                                                                                    |
|-------------------------|----------------------------------------------------------------------------------------------------|
| <b>Dilution range</b>   | The intended application should be selected according to the customer's experimental requirements. |
| <b>Specificity</b>      | No cross reactivity with other proteins.                                                           |
| <b>Cross Reactivity</b> | No cross-reactivity with other proteins.                                                           |
| <b>Antibody Type</b>    | Primary Antibody                                                                                   |
| <b>Storage</b>          | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.                   |
| <b>Note</b>             | For research use only                                                                              |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                                                    |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)