

## Product Datasheet

### CD34 Antibody (orb2623977)

|                     |                                                                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number      | orb2623977                                                                                                                              |
| Description         | CD34 Antibody                                                                                                                           |
| Reactivity          | Human                                                                                                                                   |
| Conjugation         | Unconjugated                                                                                                                            |
| Tested Applications | FC                                                                                                                                      |
| Target              | CD34                                                                                                                                    |
| Preservatives       | Purified antibody in PBS with 0.05% sodium azide.                                                                                       |
| Concentration       | Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration.                                            |
| Storage             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note                | For research use only                                                                                                                   |
| Isotype             | Mouse IgG2a                                                                                                                             |
| Clonality           | Recombinant                                                                                                                             |
| Clone Number        | 2F7                                                                                                                                     |
| Antibody Type       | Primary Antibody                                                                                                                        |
| Uniprot ID          | <b>P28906</b>                                                                                                                           |
| Entrez              | <b>947</b>                                                                                                                              |
| Dilution Range      | FCM: 1/100-1/200                                                                                                                        |
| Expiration Date     | 12 months from date of receipt.                                                                                                         |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)