

Product Datasheet

CD276(B7H3) Antibody (orb2623983)

Catalog Number	orb2623983
Description	CD276(B7H3) Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Target	CD276(B7H3)
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Recombinant
Clone Number	5D6
Antibody Type	Primary Antibody
Uniprot ID	Q5ZPR3
Entrez	80381
Dilution Range	FCM: 1/100-1/200
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)