

## Product Datasheet

### (DANRE) atg9a Antibody (orb2631152)

|                            |                                                                                                                                                                          |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | This DANRE atg9a antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 573-605 amino acids from the human region of DANRE atg9a. |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                   |
| <b>Reactivity</b>          | Zebrafish                                                                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                             |
| <b>Tested Applications</b> | WB                                                                                                                                                                       |
| <b>Immunogen</b>           | This DANRE atg9a antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 573-605 amino acids from the human region of DANRE atg9a. |
| <b>Target</b>              | (DANRE) atg9a                                                                                                                                                            |
| <b>Concentration</b>       | Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration.                                                                             |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                  |
| <b>Note</b>                | For research use only                                                                                                                                                    |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                               |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                                                         |
| <b>Uniprot ID</b>          | <b>A2RV24</b>                                                                                                                                                            |
| <b>Dilution Range</b>      | WB: 1/2000                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                          |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)