

Product Datasheet

(DANRE) atg9a Antibody (orb2631279)

Description	This DANRE atg9a antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 515-548 amino acids from the human region of DANRE atg9a.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This DANRE atg9a antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 515-548 amino acids from the human region of DANRE atg9a.
Target	(DANRE) atg9a
Concentration	Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Uniprot ID	A2RV24
Dilution Range	WB: 1/2000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)