

Product Datasheet

Penicillin Antibody (orb2637086)

Catalog Number	orb2637086
Category	Antibodies
Description	<p>This antibody reacts mainly with the thiazolidine ring of penicillin, but not with the lactam ring. May react with the fused beta lactam/thiazolidine ring of the penicilloyl group and would appear sensitive to modification of the shared N atom of the lactam/thiazolidine ring by substitution or conjugation. It appears insensitive to the structure of the side chain of the penicilloyl group. This antibody reacts with the following penicillin; Benzylpenicillin, Ampicillin, Amoxicillin and 6-Aminopenicillanic acids. Penicillin is a group of Beta-lactam antibiotics used in the treatment of bacterial infections caused by susceptible, usually Gram-positive, organisms. -lactam antibiotics work by inhibiting the formation of peptidoglycan cross-links in the bacterial cell wall, which results in cytolysis. This antibody is useful in the study of allergy to penicillin.</p>
Clonality	Monoclonal
Species/Host	Mouse
Isotype	Mouse IgG1, kappa
Conjugation	Unconjugated
Reactivity	All
Buffer/Preservatives	1 mg/ml in 1X PBS; rAlbumin free, sodium azide free
Purification	Protein A/G affinity
Immunogen	Penicilloyl-transferrin conjugate was used as the immunogen for the Penicillin antibody.
Tested applications	ELISA

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	ELISA: 1-5ug/ml (order BSA-free format for coating)
Application notes	Optimal dilution of the Penicillin antibody should be determined by the researcher.
Antibody Type	Primary Antibody
Clone Number	Pen-9
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free Penicillin antibody (clone Pen-9) as confirmation of integrity and purity.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)