

Product Datasheet

Elranatamab Antibody (orb2646426)

Catalog Number	orb2646426
Category	Antibodies
Description	Elranatamab (PF-06863135) is a humanized bispecific antibody that targets BCMA (B-cell maturation antigen) on myeloma cells and CD3 on T cells. By targeting these molecules, Elranatamab activates T cells and directs them to induce a cytotoxic T cell response against myeloma cells. This compound is used in research for relapsed or refractory multiple myeloma.
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	2408850-14-4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)