

Product Datasheet

Emirodatamab Antibody (orb2646427)

Catalog Number	orb2646427
Category	Antibodies
Description	Emirodatamab (AMG-427) is an anti-FLT3/CD3 bispecific T-cell engager with an extended half-life. It can induce high cytotoxicity in primary AML progenitor cells and enhance FLT-3 expression. Combining Emirodatamab with anti-PD-1 antibodies increases T-cell dependent cytotoxicity (TDCC). This compound is under investigation for use in relapsed or refractory AML.
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	2449199-61-3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)