

Product Datasheet

Direct Extraction Buffer (orb2646461)

Catalog Number	orb2646461
Category	Tools
Description	Direct Extraction Buffer allows an easy and fast extraction of DNA and RNA directly from blood, swabs and animal- or plant tissue. The buffer is optimized for use in combination with Direct PCR or RT-PCR master mixes like qPCR ProbesMaster (PCR-396) or SCRIPT Direct RT-qPCR ProbesMaster (PCR-528). The mix allows DNA and RNA preparation within 3-5 minutes and with a minimum of pipetting steps. It is especially recommended for:
Form/Appearance	liquid
Concentration	5 x
Storage	store at 4°C or -20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)