

Product Datasheet

RNase Inhibitor - recombinant (orb2646471)

Catalog Number	orb2646471
Category	Proteins
Description	<p>RNase inhibitor is a recombinant protein which completely inhibits a broad spectrum of eukaryotic RNases, including RNase A, B and C. It inhibits RNases by binding noncovalently in a 1:1 ratio with high affinity ($4 \times 10^{-14}M$). It does not inhibit the RNases I, T1, T2, H, U1, U2 and C13. In addition, RNase inhibitor shows no inhibition of polymerase or reverse transcriptase activity and so can be used for cDNA synthesis and in one-step RT-PCR reactions. The murine version of RNase inhibitor lacks the pair of cysteines identified in human version, therefore it has significantly improved resistance to oxidation.</p>
Form/Appearance	liquid
Concentration	40 units/ μ l
Application notes	<p>RNase inhibitor should be applied in a working concentration of 40 units (1 μl) per 50 μl reaction volume. The application of RNase inhibitor is recommended for:</p>
Storage	store at -20 °C. avoid freeze/thaw cycles
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)