

Product Datasheet

Liquid qPCR GreenMaster for Lyophilization (orb2646478)

Catalog Number	orb2646478
Category	Tools
Description	Liquid real-time PCR master mix with SybrGreen for lyophilization of qPCR assays
Concentration	2.5x conc.
Application notes	Spectroscopic Properties: dye bound to DNA: λ_{exc} 494 nm λ_{em} 521 nm
Storage	store at -20 °C. store dark, avoid freeze/thaw cycles (up to 10 freeze-thaw cycles without reduction in performance). storage at 4 °C for up to 1 week
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)