

Product Datasheet

Liquid qPCR GreenMaster for Lyophilization (orb2646479)

Catalog Number	orb2646479
Category	Tools
Description	Liquid real-time PCR master mix with SybrGreen for lyophilization of qPCR assays
Storage	store at -20 °C, store dark. storage at 4 °C for up to 1 week;. up to 10 freeze-thaw cycles without reduction in performance
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)