

Product Datasheet

AffiniPure Donkey Anti-goat IgG(H+L) secondary antibody, Biotin Conjugated (orb2646598)

Description	AffiniPure Donkey Anti-goat IgG(H+L) secondary antibody, Biotin Conjugated
Species/Host	Donkey
Reactivity	Goat
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Goat IgG (whole molecular)
Form/Appearance	Concentrated, Liquid
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Application notes	Immunohistochemistry, 5-10 µg/mL ELISA, 0.05-0.2 µg/mL
Antibody Type	Secondary Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)