

Product Datasheet

Anti-Chlamydomonas reinhardtii CrPtch1 Antibody (orb2646764)

Description	Anti-Chlamydomonas reinhardtii CrPtch1 Antibody. This antibody reacts with Chlamydomonas reinhardtii.
Species/Host	Rabbit
Reactivity	Plant
Conjugation	iFluor647
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)