


Product Datasheet

Human PP (Pepsin) ELISA Kit (orb2648021)

Catalog Number	orb2648021
Category	Assays and Kits
Description	Human PP (Pepsin) ELISA Kit
Target	PP
Reactivity	Human
Assay Type	Sandwich
Assay Time	120 minutes
Range	62.5-4000pg/ml
Sensitivity	37.5pg/ml
Sample Types	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Storage	2 - 8 °C for 12 months
Note	For research use only



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558