

## Product Datasheet

# Human Transforming Growth Factor Beta 2 (TGF- $\beta$ 2) ELISA Kit (orb2648655)

<b>Catalog Number</b>	orb2648655
<b>Category</b>	Assays and Kits
<b>Description</b>	Human Transforming Growth Factor Beta 2 (TGF- $\beta$ 2) ELISA Kit is for the quantitative determination of TGF- $\beta$ 2 levels in human serum, plasma, culture media, or any biological fluid.
<b>Storage</b>	Store all reagents at 2-8°C. Product stored under these conditions should be stable for 6 months from the production date.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)