

Product Datasheet

Rat TNF alpha Quick ELISA Kit (orb2648753)

| | |
|----------------------------|--|
| Catalog Number | orb2648753 |
| Category | Assays and Kits |
| Description | Rat TNF alpha Quick ELISA Kit |
| Target | Tumor necrosis factor |
| Reactivity | Rat |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: L80-L235 |
| UniProt ID | P16599 |
| MW | 25806 Da |
| Tested applications | ELISA |
| Specificity | Natural and recombinant rat TNF alpha |
| Cross Reactivity | There is no detectable cross-reactivity with other relevant proteins. |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Sensitivity | 1 pg/ml |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| Storage | Store at 4°C for 6 months. (Ships with gel ice, can store for up to 3 days in room temperature. Refrigerate upon receipt.) |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)