

Product Datasheet

Anti-Tomato Coatomer subunit beta Antibody (orb2648871)

Description	Anti-Tomato Coatomer subunit beta Antibody. This antibody reacts with Tomato.
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Uniprot ID	A0A3Q7IMC4
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)