

Product Datasheet

Chicken C-MYC Epitope Tag Antibody (orb2652587)

Catalog Number	orb2652587
Category	Antibodies
Description	Chicken C-MYC Epitope Tag Antibody
Species/Host	Gallus
Conjugation	Unconjugated
Form/Appearance	Provided in phosphate buffered saline with 0.09% sodium azide
Concentration	See vial for concentration.
Tested applications	ELISA, WB
Antibody Type	Primary Antibody
Source	Chickens were immunized with EQKLISEEDL (C-MYC) conjugated to KLH. After multiple immunizations in Freund's adjuvant, eggs were collected and the immunoglobulin fraction prepared. Antibodies were immunoaffinity purified using the peptide immobilized on a solid support.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)