

Product Datasheet

4-O-Caffeoyl-3-O-syringoylquinic acid (orb2653536)

Description	4-O-Caffeoyl-3-O-syringoylquinic acid
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{25}H_{26}O_{13}$
Purity	> 98% (HPLC)
MW	534.5
CAS Number	1207645-17-7
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)