

Product Datasheet

Mouse Interleukin 13,IL-13 ELISA Kit (orb2656046)

Catalog Number	orb2656046
Description	Mouse Interleukin 13,IL-13 ELISA Kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Assay Time	1-5h
Uniprot ID	P20109
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)