

Product Datasheet

Human epidermal growth factor receptor 2 (sp185/HER2) ELISA Kit (orb2656048)

Catalog Number	orb2656048
Category	Assays and Kits
Description	The Human epidermal growth factor receptor 2 (sp185/HER2) ELISA Kit detects the target v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian) in human samples. It measures levels in sample types including serum, plasma and cell culture supernates and cell lysates and tissue homogenates. The detection range is 7.8 pg/mL-500 pg/mL. The kit offers a sensitivity of 3.9 pg/mL.
Target	v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)
Reactivity	Human
UniProt ID	P04626
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	7.8 pg/mL-500 pg/mL
Sensitivity	3.9 pg/mL
Sample Types	serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)