

## Product Datasheet

# Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit (orb2656058)

<b>Catalog Number</b>	orb2656058
<b>Category</b>	Assays and Kits
<b>Description</b>	The Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit detects the target Tumor necrosis factor (TNF/TNFA/TNFSF2) in porcine samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 23.5 pg/mL-1500 pg/mL. The kit offers a sensitivity of 11.75 pg/mL.
<b>Target</b>	Tumor necrosis factor (TNF/TNFA/TNFSF2)
<b>Reactivity</b>	Porcine
<b>UniProt ID</b>	<b>P23563</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	23.5 pg/mL-1500 pg/mL
<b>Sensitivity</b>	11.75 pg/mL
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)