

Product Datasheet

FCGRT Antibody (orb2656097)

| | |
|-----------------------------|---|
| Catalog Number | orb2656097 |
| Category | Antibodies |
| Description | This FCGRT Recombinant monoclonal Antibody is an unconjugated monoclonal product. It targets FCGRT using a Recombinant Human FCGRT protein as the immunogen. This antibody is suitable for ELISA, IHC. Purification: Affinity-chromatography. |
| Target | FCGRT |
| Clonality | Monoclonal |
| Isotype | hIgG1 |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Buffer/Preservatives | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification | Affinity-chromatography |
| Immunogen | Recombinant Human FCGRT protein |
| UniProt ID | P55899 |
| Tested applications | ELISA, IHC |
| Antibody Type | Recombinant Antibody |
| Clone Number | 23H4 |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

| | |
|------------------------|---|
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)